



MANUAL BOOK AALTO

A. DENTAL MONITORING PARTNER

1.CREATE PRACTICE & DOCTOR

DentalMonitoring

Notifications

5

Patients

Doctors

To-do List

Configuration

Tools

Tammy DSO Employee (bitesoft)

Academy

Help

My Account

Logout

Practice list

Practice edition

Copy shipping info to billing info

PRACTICE NAME ("BUILDING NAME")

BILLING NAME ("COMPANY NAME")

SHIPPING E-MAIL

BILLING E-MAIL

SHIPPING STREET ADDRESS

BILLING STREET ADDRESS

SHIPPING ZIP CODE

BILLING ZIP CODE

SHIPPING CITY

BILLING CITY

SHIPPING COUNTRY

BILLING COUNTRY

SHIPPING COUNTRY REGION

BILLING COUNTRY REGION

LANGUAGE

DSO

LATITUDE

LONGITUDE

ENABLE BOOKING BY PATIENT

ENABLE PATIENT SHOP

SCANBOX *** ADAPTERS SUPPLIED

SCAN PROCESS VERSION

SCAN SEQUENCE ORDER

MANUFACTURERS

3D TOTEM

ALIWELL

ANGEL ALIGNER

ARCHFORM

CHRISTALINE

EON

GENERIC MANUFACTURER

X-LINE

LUXCREO

ONESMILE

ORTHOFX

SECRET ALIGN

SMARTEE

SMILERS

SPARK

TEST

Save

Home

Regulatory label

DentalMonitoring

Partner Monitoring System Access Manual

- 1.Click the Doctor menu.
- 2.Select Practice.
- 3.Complete the practice profile information, then click Save.

A. MONITORING PARTNER

1. CREATE PRACTICE & DOCTOR

DentalMonitoring

Notifications

5

Patients

Doctors

To-do List

Configuration

Tools

Library

DSO Employee (new)

Academy

Help

My Account

Logout

« Doctor list

Infos

New Doctor

E-MAIL

☐ DISABLED

FIRST NAME

LAST NAME

GENDER

none

ACADEMIC STATUS

none

☒ ENABLE DIRECT MESSAGING FROM PATIENTS

☐ RESEARCHER?

SECONDARY E-MAIL (ALWAYS IN CC)

OPTIONAL

LANGUAGE

Select

DSO

Select

LANGUAGE

Select

DSO

Select

EXTERNAL DOCTOR ID (INTEGRATION, PMS...)

OPTIONAL

MOBILE PHONE

FR (+33)

FEATURE SET

Advanced

☒ SEND DAILY MONITORING REPORT E-MAIL?

Doctor's main practice

Select an existing practice

Or create a new one

» Copy shipping info to billing info

PRACTICE NAME ("BUILDING NAME")

BILLING NAME ("COMPANY NAME")

CLIENT FIRST NAME

CLIENT LAST NAME

PHONE NUMBER

FR (+33)

FAX

OPTIONAL

☐ SEND PASSWORD CREATION E-MAIL

Save

Doctor Setup

1. Click the Doctor menu.
2. Select Doctor.
3. Complete the doctor profile information.
4. Select the practice that has been created.
5. Click Save.

A. MONITORING PARTNER

1. CREATE PRACTICE & DOCTOR

DentalMonitoring Notifications **5** Patients **Doctors** To-do List Configuration Tools

Tammy DSO Employee (bitesoft) Academy Help My Account Logout

« Doctor list Doctor's patient list »

Infos Timeline Agreements Login as

Arsyil Rahman [doctor]

PROFILE ID: **Z-B6A8-7FF9**

LOCAL TIME: 2026-02-07 04:59:41 (America/Aruba)

ACTIVE PATIENTS: **0** - Manage protocols - View invoice list

DM ENGAGE:

KPIS:

Practices

- Arash Klinik : +672544545 (kota jepara) - Edit

TWO-FACTOR AUTHENTICATION **Off**

Two-factor authentication provides an extra level of security when logging into your account.

E-MAIL

arsyil@gmail.com

☐ DISABLED

FIRST NAME

Arsyil

GENDER

m

Protocol Management

1. Select one doctor from the list.
2. Click Manage Protocol.
3. Click Create New Protocol.
4. Enter the protocol data, then click Save.

DentalMonitoring Notifications **5** Patients Doctors To-do List **Configuration** Tools

Tammy DSO Employee (bitesoft) Academy Help My Account Logout

Protocols Quickstarts Patient instructions Team instructions Quick replies

My activated protocols / Arsyil Rahman

Your protocols, ready for use.

API DENTIST
Last updated: 2026-01-29 15:10 GMT+7


Protocol templates

Be inspired by our selection of pre-made protocol templates, designed to address various treatment scenarios. These templates need to be reviewed and added to your protocols before they can be applied to patients.

API DENTIST
Last updated: 2025-12-10 12:28 GMT+7


A. MONITORING PARTNER

1. CREATE PRACTICE & DOCTOR





Quickstart Setup

1. Select the Quickstart menu.
2. Click New Quickstart.
3. Select a protocol and choose a treatment.
4. Click Save.



B. AALTO

1. CREATE PRACTICE




AALTO System Access


1. Log in to the AALTO system.
2. Select the Practice menu.
3. Click Create New Practice.
4. Complete the practice information.
5. Click Create Practice.

B. AALTO

1. CREATE PRACTICE



Open the Practice item, then click Assigned Practice.
Select the doctor, then click Assign.



B. AALTO

1. CREATE PRACTICE

Aalto

Dashboard

Leads7

Invoices19

Patients21

Practices

Orthodontists

Users

Messages

Settings

Search by name, email, or location...

Filter by status

Name	Address	Country	Email	Phone	Tier	Status	Created	
Jabar Hiji Klinik	jalan bandung, bandung, Jawa Barat	Cayman Islands	jabar1@gmail.com	+61324556887	bronze	Active	2/3/2026	
test practice integration	asd, asd, asd	Australia	test@practice.com	+6112345678	bronze	Active	1/30/2026	
Maman Firman	Jl. rancakalong, Sumedang, Jawa barat	Austria	coc05setiawann@gmail.com	+61323456789	gold	Inactive		
pipit practice	test, tset, tset	United States	pipit.practice@gmail.com	+187654321	platinum	Active		
eko clinic	asd, asd, asd	Austria	eko@clinic.com	+6112345678	bronze	Active		

Actions

View details

Edit practice

Deactivate

Delete practice

Practice Management

1. Select a practice.
2. Click the three-dot (:) menu.
3. Select Action.
4. Choose one of the following options:
 - View Detail to view practice information.
 - Edit Practice to update practice data, then click Update Information.
 - Deactivate to disable the practice.
 - Delete Practice to remove the practice.

B. AALTO

1. CREATE ORTHODONTIST

Orthodontists

Manage orthodontists within your provider team. Assign them to clinics, define their responsibilities, and control their access to clinic cases.

Search by name, email, or location... Filter by status

	↑↓ Name	Phone	↑↓ Status	↑↓ Created	
SJ	Sarah Johnson sarah.johnson@ortho.com	+62811222555	Active	1/13/2026	⋮
LA	Lisa Anderson dr.lisa@ortho.com	+62811222666	Active	1/13/2026	⋮
IO	Ikhsan Orthodontist ikhsan@orthodontist.com	0897654321	Active	1/13/2026	⋮

Actions

- View details
- Deactivate
- Delete doctor

1. Select the Orthodontist menu.
2. Click Create New Orthodontist.
3. Complete the orthodontist information, then click Create Orthodontist.
4. Click the three-dot (⋮) action menu:
5. View Detail to display orthodontist information.
6. Deactivate to disable the doctor.
7. Delete Doctor to remove the doctor data.

Create New Doctor

Add a new orthodontist to the system

Personal Information

First Name * Last Name *

John Doe

Phone Number

+61 +61 Phone Number

Account Credentials

Email Address *

orthodontist@example.com

Password * Confirm Password *


Enter a password. Must contain:

- ⓧ Password must be at least 12 characters
- ⓧ Password must contain at least one number
- ⓧ Password must contain at least one lowercase letter
- ⓧ Password must contain at least one uppercase letter
- ⓧ Password must contain at least one special character

Cancel Create Orthodontist


B. AALTO

1. CREATE USER



User Actions

1. Select a user from the list.
2. Click the three-dot (:) action menu.
3. Choose one of the following options:
4. View Details to display user information.
5. Edit User, update the data, then click Update Information.
6. Deactivate to disable the user.
7. Reset Password to change the password. The reset link will be sent to the user's registered email address.
8. Delete User to permanently remove the user data.



User Management

1. Select the User menu.
2. Click Create New User.
3. Complete the user information.
4. Select the role and access level.
5. Select a doctor from Dental Monitoring.
6. Click Create User.

B. AALTO

1.CREATE LEAD PATIENT (ACCOUNT ADMIN)

Leads Management

Manage and track patient leads assigned to your practice.

Assigned

Rejected

Accepted

Search by name, email, or location...

Assigned

Name	Status	Email	Phone Number	Appointment Date	Assigned Practice
test create lead ikhsan	Waiting Confirmation	j7vyh@virgilian.com	+61412345678	N/A	—
eko haryadi	Waiting Confirmation	ekoharyadi416@gmail.com	+6285320080416	N/A	Echo Practice bronze
eko haryadi	Waiting Confirmation	ekoharyadi416@gmail.com	+6285320080416	N/A	Echo Practice bronze

View Details

Cancel Lead

- 1.Select the Lead Patient menu.
- 2.Click Create Lead.
- 3.Complete the lead information and medical history.
- 4.Click Create Lead to save.

Create New Lead

Practice Assignment

Practice

Assigned Practice

Personal Information

First Name

Last Name

Email Address

Phone Number

Address Information

Street Address

City

State/Province

Zip/Postal Code

Country

Medical History

Please indicate if you currently have or have ever had any of the following conditions

Heart condition Diabetes Autoimmune condition Bleeding disorder

Exposure to asbestos Osteoporosis Blood cancer TB or jaw issues

Chronic pain disorder None of the above

Other medical conditions (please specify)

Please provide additional details

Are you currently taking any medications (including supplements)?

No Yes - please list

Please provide additional details

Do you have any known allergies (including latex or dental materials)?

No Yes - please list

Please provide additional details

Are you currently pregnant or planning pregnancy?

No Yes Not applicable

Dental & Orthodontic History

Please provide information

Have you had orthodontic treatment before?

No Braces Clear aligners

If yes, please provide details (type, location, duration)

Please provide additional details

Which of the following best describes your current concern?

Oral health Bleeding / pain Swelling Underbite Crooked bite

Loose / shifting or clicking Teeth grinding or wear None of the above

Have you had any of the following dental treatments?

Crown Veneer Dental implants Root canal treatment Endodontic therapy

None of the above

Lead Details

Practice Assigned (this will be used for reporting)


Additional Notes

Please provide additional comments or special requests

Cancel Create Lead

B. AALTO

1. CREATE PATIENT (ACCOUNT DOCTOR)



The 'New Patient' form is divided into two main sections: 'Basic Information' and 'Patient Information'. The 'Basic Information' section includes fields for 'First Name*', 'Last Name*', 'Date of Birth*', and 'Gender*'. The 'Patient Information' section includes fields for 'Email Address*' and 'Phone Number*'. There are also sections for 'Dentist Assignment' and 'Diagnostic materials'. A black arrow points from the '+ Create New' button in the dashboard to the 'NEW CASE' button in the form.

New Patient

Basic Information
Enter the patient's personal and contact information to create their case.

Dentist Assignment
Assign a Dentist to this case if not already assigned.

Select Dentist*

Search and select a dentist...

Patient Information
Personal details and demographics

First Name*
Enter first name

Last Name*
Enter last name

Date of Birth*
dd-....yyyy

Gender*
Male

Contact Information
Email and phone details

Email Address*
john@example.com


Phone Number*
+1 Phone Number

Address Information
Shipping and mailing address

1. Create Patient Case (Doctor Account)
2. Select the Patient menu.
3. Click Create New Patient.
4. Complete the following information:
5. Basic Information
6. Medical History
7. Treatment Goal
8. Diagnostic Materials
9. Click Create Patient.

B. AALTO

1.ACCEPTED LEAD & BOOK APPOINTMENT



- 1.Accept Lead Patient (Doctor)
- 2.Log in as Doctor.
- 3.Select Lead Patient.
- 4.Click Accept Lead.
- 5.Confirm acceptance.
- 6.Select the doctor and booking date.
- 7.Click Booking Appointment.
- 8.Complete the Treatment Goal and Diagnostic Materials.
- 9.Click Update Information.

B. AALTO

1. APPROVAL BY ORTHODONTIST

The screenshot shows the Aalto Case View interface. The main content area displays the 'Request Clarification' modal, which is a form for providing details on what needs to be clarified. The modal has a 'Notes for the dentist' section with a text input field and a 'Send Request' button. The 'Need Clarification' button is highlighted with a red arrow pointing to the 'Notes for the dentist' section. The 'Write Summary' button is also highlighted with a red arrow pointing to the 'Write Summary' modal on the right.

Aalto CASE VIEW

CASE BSC 20260207 00002 - FEB 07, 2026
Jaja Miharja

Treatment Plan Summary

Patient Information

Basic info Diagnostic materials Treatment goals

Case Detail
View the detailed information for this case.

Dentist Assignment
Assign a dentist to this case if not already assigned.

Select Dentist *

Search and select a dentist...

Patient Information
Personal details and demographics

First Name
Enter first name

Date of Birth
dd-yyy

Contact Information
Email and phone details

Email Address
jaja@example.com

Phone Number
+1

Request Clarification

Notes for the dentist
Please provide details on what needs to be clarified.

Cancel Send Request

The screenshot shows the 'Write Summary' modal form. It contains sections for uploading a treatment plan PDF, selecting a case type, and providing an overview of the treatment plan. The 'Write Summary' button is highlighted with a red arrow pointing to the 'Write Summary' modal on the right.

Write Summary

Upload treatment plan (PDF)
Please upload the treatment plan PDF in the field below

Drop items here or Browse files

Select Case Type *

Express

Treatment plan overview
Please provide an overview of the treatment plan for the dentist by filling out the forms below.

What you can achieve with this treatment

Watch out for

Why did I choose this approach?

Inform your patients

Cancel Submit Summary

1. Select a case to open the Case View.
2. If clarification is required, click Need Clarification and send a request.
3. Click Approve Treatment or Reject.
4. If all data is complete and correct, select Write Summary.

B. AALTO

1. BOOK APPOITMENT

The screenshot shows the 'CASE VIEW' interface for a patient named Jaja Miharja. The patient's profile includes their date of birth (1998-02-04), email (jajamiharja@gmail.com), phone number (+62 873538382), and home address (jalan jalan men). The case status is 'Submitted Case', with a 'Submitted' button and a 'View Submission' link. The 'Case tracker' section shows a 'Case submitted' status with a note: 'Case has been submitted for review • February 7, 2026 at 10:53 PM'. The 'Treatment Plan Summary' section is currently empty.

The screenshot shows the 'EDIT CASE' interface. The 'Edit Case' section includes three expandable panels: 'Basic info' (Jaja Miharja, (28 y/o)), 'Treatment goals' (Choose which practical features you'd like to use for this case...), and 'Diagnostic materials' (Our team review these diagnostic materials and check for any issues that could complicate treatment). The 'Treatment Goals' section includes three questions: 'Do you prefer to use IPR? Standard test *' (Yes, No, No, except to improve black triangle), 'Do you prefer to use attachments? *' (Yes, No), and 'Do you prefer to use bite ramps? *' (Yes, No).

1. Resubmit Case Information (Doctor)
2. Log in to the Doctor account.
3. Open the case that requires resubmission.
4. Complete the requested information.
5. Click Update Information.

B. AALTO

1. APPROVE BY ORTHODONTIST

CASE VIEW

CASE ESC-20260207-00002 - FEB 07, 2026
Jaija Mharja

test dentist practice

Treatment Plan Summary

Patient Name	Case Number	Case Type	Created by	Status
Jaija Mharja	ESC-20260207-00002	Misc	to: test orthodontist practice	Draft

Version
Version 1 (Draft) at Feb 8, 2026, 10:42 AM
Created: Feb 08, 2026, 10:42 AM
Updated: Feb 08, 2026, 10:42 AM

Treatment Plan Document
(7375522-487958.pdf)

Treatment Plan Overview
What you can achieve with this treatment
ofteoff
Watch out for
tglels
Why did I choose this approach?
glafld dgtf
Inform your patients
tfgd

- REC's
- Staging
- Telreport
- Attach
- Retic
- Super
- Sild
- Roots
- Num
- PRE
- Space
- Oc P
- TMT
- MDA
- Button
- OUOB
- Arch
- CSHME
- Mess
- Lang

Name: 咬合正
ID: T11142
Plan: Stage Plan 1A-1A

Initial Phase
Stage Plan 1A-1A Approved

Review Approved

This treatment plan is still in draft and has not been sent to the dentist.

Edit Plan Send to Dentist

Patient Information

Basic info Diagnostic materials Treatment goals

Case Detail
View the detailed information for this case.

Dentist Assignment
Assign a dentist to this case. If not already assigned.

Select Dentist *
test dentist

Patient Information
Personal details and demographics

First Name Last Name
Jaija Mharja

Date of Birth Gender
04-Feb-1998 male

Contact Information
Email and phone details

Email Address Phone Number
jaijamharjad@gmail.com +62 873538382

Address Information
Patient home address

Street Address City State/Province
jalan pan-mendali gunung kota city provincial pr

Appt/Suite Country
3101 SA Germany

Zip/Postal Code Relationship
34331 e.g., Spouse, Parent


Emergency Contact
Secondary contact information

Contact Name
Enter contact name

Phone Number
Relationship

1. Review by Orthodontist
2. Log in as Orthodontist.
3. Review the updated case.
4. If approved, click Send to Dentist.
5. Case status can be tracked in the Case Tracker as Treatment Plan Approved and Review Approval.




1. APPROVE BY ORTHODONTIST



- Case status can be tracked in the Case Tracker as Treatment Plan Approved and Review Approval.

B. AALTO



1. AALTO AI SUMMARIES



1. Generate Summaries Using AI
2. Treatment Summaries
3. Select Treatment Summaries.
4. Complete the required information.
5. Click Generate.
6. The AI-generated treatment summary will appear.
7. Click Save to Profile or Copy.

B. AALTO



1. AALTO AI SUMMARIES



1. Insurance Summaries
2. Select Insurance Summaries.
3. Complete the required information.
4. Click Generate.
5. The AI-generated insurance summary will appear.
6. Click Save to Profile or Copy.

B. AALTO


1. AALTO AI SUMMARIES




- PMS Summaries
- Select PMS Summaries.
- Complete the required information.
- Click Generate.
- The AI-generated PMS summary will appear.
- Click Save to Profile or Copy.

B. AALTO

1. SIGN PATIENT CONSENT DOCUMENT




Patient profile card for Jaja Miharja. The card displays personal information: Date of Birth (1998-02-04), Email (jajamiharja@gmail.com), Phone Number (+62 873538382), and Home Address (jalan jalan mendaki gunung). Below the profile information, there are two sections: 'Consent Patient' and 'Medical History'. The 'Consent Patient' section shows a 'Sent' status and a 'Patient Consent' button, which is highlighted with a red box. The 'Medical History' section shows a 'Sent' status and a 'Medical History' button, which is also highlighted with a red box.



Patient Agreement document titled 'Aalto Clear Aligner Treatment Risks & Informed Consent'. The document includes an 'Important Notice to Treating Dentist' and a 'SIGN NOW' button.



Patient Signature screen. It displays the patient's name (Jaja Miharja), the date (02/08/2026), and the patient's signature. The screen includes a 'REDRAW' button and a 'NEXT' button.




Document signed confirmation screen. It displays the date (02/08/2026) and a confirmation message: 'Document has been signed!'. The screen includes buttons for 'SEND COPY VIA EMAIL' and 'DOWNLOAD', and a footer indicating it is powered by DocuSeal - open source documents software.

1. Click the three-dot action menu.
2. Select Patient Consent.
3. Scroll through the document.
4. Click Sign In.
5. Draw the signature and click Next.
6. Set the date and click Complete.
7. The document is successfully signed and can be emailed or downloaded.

B. AALTO

1. SIGN MEDICAL HISTORY AGREEMENT



Jaja Miharja

Date of Birth
1998-02-04

Email
jajamiharja@gmail.com

Phone Number
+62 873538382

Home Address
jalan jalan mendaki gunung

Consent Patient


Sent Medical History

Awaiting patient's consent to Patient Consent

Medical History

Sent

Awaiting patient's medical history to proceed.



DocuSeal

Medical History Agreement V2


DECLINE DOWNLOAD

Aalto

Preferred

Type here...

NEXT



Section 2 — Medical History

Please indicate if you currently have, or have ever had, any of the following conditions
(tick all that apply):

☐ Heart condition

☐ Diabetes


☐ Autoimmune condition

☐ Bleeding disorder

...

NEXT

Other medical conditions (please specify):



Date

02/08/2026

Document has been signed!

SEND COPY VIA EMAIL


DOWNLOAD

Powered by DocuSeal - open source documents software

1. Click the three-dot action menu.
2. Select Medical History.
3. Scroll and review the medical history document.
4. Click Sign In.
5. Draw the signature and click Next.
6. Set the date and click Complete.
7. The document is successfully signed and can be emailed or downloaded.

B. MANUFACTURE

1.ASSIGN LAB & MANUFACTURING STATUS



1. Log in to the Manufacturing System.
2. Select the Case menu.
3. Choose a case.
4. Click Assign Lab and select a lab.
5. Click Assign Lab to confirm.
6. Update the status from Received to In Production, then click Update Status.
7. Click Save.
8. Update the status again to Shipped.
9. The case status can be tracked until it reaches Shipped.

B. MANUFACTURE

1. CASE TRACKER

Case tracker

Case submitted

Case has been submitted for review • February 7, 2026 at 10:53 PM

Treatment plan in progress

A licensed orthodontist is currently building the treatment plan. Plans typically take between 24-48 hours. • February 8, 2026 at 10:20 AM

Materials rejected

hello plese complete • February 8, 2026 at 10:21 AM

Treatment plan ready for review

Your treatment plan is ready for review. Approve or request changes. • February 8, 2026 at 10:43 AM

Treatment plan approved

Case sent to manufacture • February 8, 2026 at 12:33 PM

Production started

Case sent to manufacture • February 8, 2026 at 12:33 PM

Aligner kit in production

Production can take around 10 business days to be manufactured and shipped out. • February 8, 2026 at 04:09 PM

Central Diagnostic Lab

Quality control completed

Quality control check has been completed successfully. • February 8, 2026 at 04:09 PM

Central Diagnostic Lab

Aligner kit shipped

Your aligner kit has been dispatched and should be with you within 5 days. • February 8, 2026 at 04:10 PM

Central Diagnostic Lab



test orthodontist practice

Orthodontist

- The case status can be tracked until it reaches Shipped.

B. MANUFACTURE


1. DELIVERY & TREATMENT STATUS



1. Open the case.
2. Scroll to the bottom.
3. Click Mark as Delivered and confirm.
4. Status will change to Start Treatment.
5. Click Start Treatment.
6. Once finished, click Completed Treatment.

B. MANUFACTURE

1. START TREATMENT



- 1. Open the case.
- 2. Scroll to the bottom.
- 3. Click Mark as Delivered and confirm.
- 4. Status will change to Start Treatment.
- 5. Click Start Treatment.
- 6. Once finished, click Completed Treatment.

B. MANUFACTURE

1. OPEN MONITORING



A patient profile card for Jaja Miharja. It features a circular profile picture of a woman with blonde hair. Below the photo is the name "Jaja Miharja". Further down, it lists personal details: "Date of Birth 1998-02-04", "Email jajamiharja@gmail.com", "Phone Number +62 873538382", and "Home Address jalan jalan mendaki gunung". At the bottom, there is a rounded button labeled "Open Monitoring".

Jaja Miharja


Date of Birth
1998-02-04

Email
jajamiharja@gmail.com

Phone Number
+62 873538382

Home Address
jalan jalan mendaki gunung

Open Monitoring



A screenshot of a "Patient Monitoring" setup screen for Jaja Miharja. The screen is divided into several sections. On the left, there's a "Monitoring Information" sidebar with a "Start Monitoring" button. The main area is titled "Patient Information" and contains various input fields for personal and medical data. At the bottom, there are checkboxes for notification preferences and a "Doctor" field.

Patient Monitoring
Jaja Miharja

Monitoring Information

Plan
Started
Patient App: Not Activated
Scan Box
Next Scan
Frequency
Upper/Lower
Aligner #
Excluded Teeth

Start Monitoring

Notifications to Review

There is no notification to review

Patient Information
Basic information and settings

First Name: Jaja
Last Name: Miharja
Email: jajamiharja@gmail.com
Secondary Email (Optional):
Language: English
Date of birth: 04-Feb-1998
Push notification alert: 19:00
Mobile phone (Optional): +62 873538382
Practice Software Id (Optional):
Scan Process Version: Scan process V2
Scan Sequence Order: Not specified, Aligners in first, Aligners in last
Enable Push Notification: ☒ ☐ Enable Email Notification
Adapter Instruction: ☐ Has Led Adapter: ☐


Doctor
Practice

1. Monitoring
2. Open the case.
3. On the left profile menu, click Open Monitoring.
4. Monitoring information will be displayed, including: Information, Agreement, Public Profile


6. Click Start Treatment to begin monitoring.
7. Select Quickstart Monitoring.
8. Choose Pause or Stop Monitoring if needed.
9. Click Completed Treatment when finished.

B. MANUFACTURE

1.OPEN MONITORING




1. Click Start Treatment to begin monitoring.
2. Select Quickstart Monitoring.
3. Choose Pause or Stop Monitoring if needed.
4. Click Completed Treatment when finished.



B. INVOICE

1. INVOICE



Create Invoice

Invoice Information

Currency *

Australian Dollar (A\$)

Due Date *

Feb 8, 2026

Bill To

Practice Name

test practice integration

Email

test@practice.com

Phone Number

+6112345678

Address

asd

Items Invoice

Root Canal Treatment

Search...

Root Canal Treatment

\$300.00

Dental Crown

\$500.00

Teeth Whitening

\$200.00

No Items Added

There are no invoice items added yet. Please add items to proceed.

Cancel

Send Invoice

Add Custom Invoice Item

Item Name *

Enter item name

Quantity *

1

Price *

0


Cancel

Add Item

- Create Invoice (Admin)
- Log in as Admin.
- Select Invoice.
- Choose a case in Invoice Management.
- Complete invoice details and add items.
- Click Create Invoice.
- The invoice will appear in the case.

B. INVOICE

1. INVOICE



Invoice Preview

Invoice Information

Currency*

Australian Dollar (A\$)

Due Date*

Feb 8, 2026

Bill To

Practice Name

test practice integration

Email

test@practice.com

Phone Number

+6112345678

Address

asd

Items Invoice

Item Name	Quantity	Price	Total
Root Canal Treatment	1	\$300.00	\$300.00
Total			\$300.00

Cancel

Pay

- 1.Pay Invoice (Doctor)
- 2.Log in as Doctor.
- 3.Select Invoice.
- 4.Choose a case.
- 5.Click Action → Pay.
- 6.Review the invoice preview and click Pay.

B. INVOICE

1. INVOICE

Aalto sandbox

Sandbox

A\$300.00

Due February 8, 2026 Due today

↓ Invoice PDF

To

test practice integration

From

Aalto sandbox

Invoice

#DXYY5PUD-0002

Memo

Invoice for Case Number : BSC-20260207-00002

View invoice and payment details >

Choose currency for payment

🇮🇩

IDR 3,690,331.95

🇺🇸

A\$300.00

1 AUD = 12,301.1065 IDR

Choose how you'd like to pay

📁

Saved

VISA

Visa

**** 4242

⋮

🗳️

Card

Notwithstanding the logo displayed above, when paying with a co-branded eftpos debit card, your payment may be processed through either card network.

Pay

Powered by stripe

Terms Privacy

Invoice paid

A\$300.00

View invoice and payment details >

Charge amount

IDR 3,690,331.95

Invoice number

DXYY5PUD-0002

Payment date

February 8, 2026

Payment method

Visa **** 4242


Download invoice

Download receipt

- Select Card as the payment method.
- Enter card details and click Pay.
- After successful payment, payment details will be displayed and can be downloaded.

B. MESSAGE

1. MESSAGE



- Messaging
- Log in as Doctor.
- Open a case.
- Click Live Message to communicate with the Orthodontist or Admin.